



O T O P L A S T I K + L A B O R T E C H N I K






## Hearing Protection



### Fixing Clip for Rubber Cord

black	864944
red	864947
blue	864948



## Care Products



### **Cedis Cleansing and Drying Set for Slim Tubes and Open Fit Hearing Aids**

Open fit hearing aids require special cleansing and drying products. To this end, *Cedis* also offers the best solution:

The comprehensive care set, packed in an attractive transparent bag, for the cleaning and drying of open fit hearing aids. Suitable for all kinds of established slim tubes!

#### **Standard**

87416

With the OtoFloss special cleansing threads your slim tubes will be dried effectively. Before this step, the cleansing tablets with active oxygen will provide an optimum cleaning. Drying capsules (with color indicator) and the drying container for the drying of hearing aids complete the set.

*Contents:* 8 cleansing tablets, 1 cleansing container with sieve, 10 OtoFloss mini special cleansing threads for slim tubes with  $\varnothing \leq 0.9$  mm, 4 drying capsules, 1 drying container

Also available as premium version!



#### **Premium**

87417

The *Cedis* Care Set Premium with Dry Case\*.

In addition to the *Cedis* cleansing products, this set also contains the electrical dry case for the comfortable drying and storing of BTE hearing aids with slim tubes.

*Contents:* 1 *Cedis* dry case with power supply, 8 cleansing tablets, 1 cleansing container with sieve, 100 OtoFloss special cleansing threads for slim tubes with  $\varnothing \leq 0.9$  mm



\* Dry case in new design (for more information go to next page) will be available from the first quarter of 2010.

**Cedis Dry Case for Hearing Systems, 12 V**

87470

Drying and storing all in one! Electrical dry case (empty) – **now with time control and automatic switch-off in power supply!** For the fast and comfortable drying and storing of hearing aids and earmolds. Regular use contributes substantially to expanding the life expectancy of hearing aids – we recommend 3 hours per use. Maximum temperature of 50 °C (122 °F). Low electricity consumption of 1 W. With light-emitting diode during the drying process. Dry case is operated with a power supply unit (12 V) which is included. Can be used with mains voltage from 100 V up to 240 V!

Available from the first quarter of 2010.



**Cedis Drying Cushion ETK 3, 12 V**

87485

Electrical drying cushion (empty) for the fast and convenient drying of all types of hearing aids and earmolds. **Now with time control and automatic switch-off in power supply.** The built-in heated loop ensures a constant and effective drying temperature. Regular use contributes substantially to expanding the life expectancy of hearing aids – we recommend 3 hours per use. Incl. temperature control. The cloth cover is removable and washable at 30 °C (86 °F). Can be used with mains voltage from 100 V up to 240 V!


Available from the first quarter of 2010.



**Cedis med WB, water-based**

The unique hearing system care for hearing aids and earmolds with long-term antimicrobial effectiveness! **Now water-based without alcohol.** The disinfectant solution protects against microbial deposits for 5 days.

- highly effective
- long-term effective
- antibacterial



Apply once, enjoy 5 days of sanitary certainty

**Cleansing Stick, 15 ml**

86840

**Cleansing Wipes, 15 pieces**

83844



**Cedis Cleansing Container, 150 ml****80200**

For cleaning BTE earmolds with *Cedis* cleansing tablets. **With integrated withdrawal device** for sieve!

**Container**, incl. sieve and lid

**Cedis Cleansing Concentrate, 250 ml****87155**

Cleans reliably and removes persistent cerumen deposits from BTE earmolds. Cleansing concentrate in a practical bottle with integrated dosing help! To be used in *Cedis* cleansing container (REF 80200). To be mixed with water to a ratio of 1:4. Produces 1.25 l cleansing solution.

**Cedis OtoFloss, mini****78039**

Special cleansing threads for sound tubes with different inner diameters, ideal for open fit hearing aids. **Now available in small package.**

**OtoFloss mini** for inner diameter  $\leq 0,9$  mm

**Bag** with 10 pieces

**Cedis MultiTool****78041**

Special tool with magnet for easy removal of batteries. The integrated brush and synthetic thread are ideal for removing persistent deposits in slim tubes.

Transparent pillow pack with 5 MultiTools



**Cedis Profi UVC Clean**

86256

The professional cleansing device – guaranteed to quickly and thoroughly sanitize the entire hearing system in your acoustician shop.

In the *Cedis Profi UVC Clean*, the surface of the hearing aids and earmolds is treated with ultraviolet light in a dry, natural and uncomplicated procedure. Viruses and bacteria are effectively destroyed through this technology, which is also used to sanitize drinking water. The high-gloss, metallized chamber is evenly illuminated due to the ring-shaped UVC lamp.

- Easy access to lamp chamber due to hinged lid.
- 3- or 5-minute timer control with security switch-off when lid is opened.
- Blue LED display indicates operation time, that is the finishing of the process.

The device is operated with a 12 V power supply unit (included). Usable for all mains voltages from 110 to 240 V.

**Cedis [my.name]**

Your name as a brand – with your own care range

Attractive – Individual – Flexible

**Cedis [my.name] Eco:**

The favorably priced starting program

Minimum order quantity: **Now 50 pieces**

Ideal for testing new products

No contractual commitments

No additional charges, but with sliding-scale prices

Attractive design

Packaging with basic information (contents, CE, etc.) in four languages:

German, English, French, Spanish

Including satisfaction guarantee

And this is how it works:

- Select products
- Select design
- Define color
- Send your logo to us

... and there it is: your professional care program with your own name.

For your individual order form, please contact us under: +49 831 5811360



**Cedis [my.name] Standard:**

For customers with an existing care program

You are already successfully selling care products on the market.

You require individual products, adjusted to your company appearance.

You want high quality with flexibility.

Our offer: **Cedis [my.name] Standard**

- Custom-made boxes **now from 250 pieces**
- No contractual commitments
- Stickers on neutral boxes with personalized design

Please ask for your individual offer under: +49 831 5811360

**Cedis [my.name] Premium:**

(Nearly) everything is possible.

You need care products in large numbers.

You want a personalized design, adjusted to existing products.

You wish to offer a different box size.

You intend to amplify your product range with products besides standard.

Our offer: **Cedis [my.name] Premium**

- **Now available from 1000 pieces**
- Boxes are printed directly in all desired languages
- Free storing of products and packaging
- We will help you with the design

For detailed information, please contact us under: +49 831 5811360



# Impression Technique



## Impression Tweezers

45630

With round tips for safe positioning of impression pads, with sliding clasp lock mechanism.



# Service and Repairs



## The Shop on Wheels

Home Visit. Hearing protection field work. Professional education.

Plastic bag?

We offer you something better!

The chic otoplastastic trolley. Your mobile packing station for everything having to do with earmolds. We will be pleased to supply you with the completely equipped trolley.



**Otoplastic Trolley (empty), aluminum black****66005****Otoplastic Trolley "Standard"****66006**

Offers many assembling possibilities! In order to ensure you have everything you might need, we chose the contents according to the high standard needs of an acoustician's master exam in Germany:

Contents:

Heine Mini 3000 Otoscope 2.5 V halogen, black	Black Foil, 10 pieces, 4.5 cm
Heine Ear Light with Acrylic Tip, black, with the new Heine Mini C clip lamp	LP-Cleaner, 250 ml
Tube Widening Pliers	Silicone Vent Tube, 0.8 mm
Scissors, round tip	Silicone Vent Tube, 1.5 mm
Bayonet Tweezers	Silicone Vent Tube, 2.0 mm
egger A/soft (8 double cartridges)	Tweezers, fine, curved
Inject Pistol	Drill Set for hard earmold materials
Mixing Canulas, 25 pieces	Bur, round, 3.5 mm
Clear Former – Impression Film (2 x 10 pieces)	Carbide Face-Plate-Trimmer
Cotton and Sewing Thread for manufacturing impression pads	Soft Polishing Disc, green, not assembled (10 pieces)
Cedis Spray with atomizer, 100 ml	Mandrel for Soft Polishing Disc
Bottom Part (for gasket forms), 40 mm	Corundum Paper, grit 180
Gasket forms, 40 mm, 2 pieces	Uniglu, 10 g with dosing tip, small, 0.5 mm
Duplicating Silicone transparent "Express", 8 cartridges	Universal Knife
egger Acryl/B Monomer (liquid), 250 ml	Spare Blade for Universal Knife (1 piece)
egger Acryl/B Polymer (powder), transparent, 0.5 kg	Scalpel with 1 Blade, flat handle
Mixing Cup, 70 ml	Spare Blades for scalpel with flat handle (5 pieces)
Mixing Spatula, small	Cotton Swabs, 200 pieces
Immersion Skewer for waxing of 1 impression (only in standard version)	Safety Visor
Dipping Wax, green, 167 g (only in standard version)	Face Mask, double-fleece, 5 pieces
Sound Tube PVC, 2.0 x 3.1 mm, transparent	Stetoclip, black
egger LP/H beige, 150 g tube	Mixing Slab, plastic
egger LP/H Modeling Paste beige, 20 g	Drill Stand Attachment with old-away lid
Dipping Jar, 60 ml	1 Set of Training Earmolds including sound tubes (the set includes each 5 x right / 5 x left)
	Paper Tissue Dispenser, 150 tissues
	Digital Timer
	Edding, fine, black

\* Laboratory devices are not included, as they are generally provided by the schools while in training or during exams.

**Otoplastic Trolley "Expert"****66007**

This trolley, in addition to the products of the "Standard" version, is equipped with machines for the easy manufacturing of transparent earmolds of UV materials:

egger EL 3, 230 V  
Dipping Wax, green, 500 g  
Wax Immersion and Duplicating Device, 230 V  
Dipping Sieve  
Digital Thermometer  
Aquapress



**Training Earmolds** including sound tubes

(the set includes each 5 x right /5 x left)

The shape of the training earmolds is based on the main characteristics of more than ¾ of all earmolds produced daily. A great auxiliary for practicing fittings and for exams! Pre-bent sound tubes 1.5 x 2.5 mm are inserted!

**19930**Composition:

- 1 earmold without venting
  - 1 earmold with venting of 1.0 mm
  - 1 earmold with venting of 1.4 mm
  - 1 earmold with venting of 2.0 mm
  - 1 earmold with venting of 2.5 mm
- (left and right respectively)



We also offer training courses on earmolds and hearing protection for technicians in our fully equipped training laboratory!

For further information, please consult our product catalog page 202.

**Stetoclip Kit Premium**

The premium composition in a single kit.

Stetoclip with elastic earplugs of soft plastic, including listening tube with top. To avoid mix-ups, every stetoclip is equipped with a name tag. Color: blue.

**97030****Earmold Vacuum Cleaner**, manually operated

Ideal for the professional removal of cerumen by creating a vacuum. Ergonomically shaped with operable release button. The tip is equipped with an exchangeable listening top (REF 88121) that can be purchased separately.

**97070**

# Presentation / Promotion Assistance

## Deco Column in Metallic Design

Light reflection draws an interesting play of colors. Delivery includes white deco ear.

<u>Color</u>		with left ear	with right ear
<b>silver-metallic</b>	small	<b>87170-L</b>	<b>87170-R</b>
	medium	<b>87171-L</b>	<b>87171-R</b>
	large	<b>87172-L</b>	<b>87172-R</b>
<b>mango-metallic</b>	small	<b>87173-L</b>	<b>87173-R</b>
	medium	<b>87174-L</b>	<b>87174-R</b>
	large	<b>87175-L</b>	<b>87175-R</b>



## Deco Display Colorline

With Colorline you can add even more color to the presentation of hearing systems. Choose from the standard colors red, black or white. The displays are also available in many other colors with a minimum order quantity of 50 pieces. Colorline is made of plastic and is available in two sizes to be able to hold either one or two ears.

You have special requests? – We can also offer individualized solutions. Please order the matching egger deco ears separately.

Dimensions:

W 18 x H 15 cm

	for 1 deco ear	for 2 deco ears
<b>black</b>	<b>88300</b>	<b>88301</b>
<b>white</b>	<b>88302</b>	<b>88303</b>
<b>red</b>	<b>88304</b>	<b>88305</b>



### Magnetic Glass Board

A multi-functional item among the presentation products with several possible uses! Experience unlimited creativity and flexibility in your product presentation set-up. Whether you want to simply attach an egger Magic Ear to present a new hearing system or make the *Cedis* care products the center of attention: You can position many products rapidly and easily by using magnets. In addition, you can label the magnetic glass board with dry-erase markers. Our magnetic glass boards combine simple design with color effects and light effects due to reflecting glass. They are available in trendy and timeless colors. You can mount the magnetic boards horizontally or vertically to the wall or set them up vertically anywhere you like (e. g. in a display window). Please order the matching egger Magic Ears separately.

Delivery includes:

Mounting set for wall  
5 magnets

Color / size (W x H)	50 x 30 cm	82 x 12 cm	80 x 50 cm
<b>green</b>	<b>83850</b>	<b>83851</b>	---
<b>red</b>	<b>83852</b>	<b>83853</b>	---
<b>white</b>	<b>83854</b>	<b>83855</b>	---
<b>black</b>	<b>83856</b>	<b>83857</b>	<b>83858</b>

Accessories:

#### Sets with 5 magnets

Special magnets with an extra high holding power for use with the magnetic glass boards.

#### Ball-Shaped, Ø 12.7 mm

- Holding power on the magnetic glass boards: 5 sheets DIN A4
- Holding power on iron: approx. 3.2 kg!

83859



#### Cubic

The cube has 2 strong sides (magnetic poles) and 4 weaker sides.

- Holding power on the magnetic glass board (strong side): 6 sheets DIN A4
- Holding power on the magnetic glass board (weaker side): 3 sheets DIN A4
- Holding power on iron (strong side): approx. 4.1 kg!

838591



Dimensions: 10 x 10 x 10 mm

#### Cylindrical

- Holding power on the magnetic glass board: 5 sheets DIN A4
- Holding power on iron: approx. 3.9 kg!

838592



Dimensions: 10 x 10 mm

#### Support Bracket for free set-up

For the magnetic glass board

**50 x 30 cm**  
**82 x 12 cm**

**83870**  
**83871**

#### Dry-Erase Markers

They always provide the right color for a colorful labeling of the egger magnetic glass board. Dry-erase markers can be wiped off without water.

**neon-blue**  
**neon-green**  
**black**  
**white**

**83860**  
**83861**  
**83862**  
**83863**



**Presentation Tray, black**

High-quality hearing systems deserve a professional presentation during a sales talk. The surface is made of synthetic leather. Pictured presentation products are not included.

Dimensions:

W 27 x D 17 cm

**87890****Bulk Bin / Sales Counter Display**

Ideal for placing give-aways on the sales counter. Stylish sales counter display in glass-clear and high-gloss design. Individual labeling with signs or labels on request. The display is delivered in two pieces that can be mounted in no time.

**87986****Poster Clamp Frame**

The matching frame for a nice poster. The offered poster clamp frames are ideal for presenting posters hanging in a display window, on the wall or from the ceiling in the sales area. The frame, with its aluminum design, is easily mounted. Posters are quickly exchanged. They are held in place by a clamp mechanism. Suitable for DIN A1 posters (841 x 594 mm).

**83890**Outer dimensions of the frame:

H 922 x W 675 mm

**Poster Clamp Strip**

Ideal when there is no time! The simple and stylish solution for attaching and hanging up posters up to 600 mm wide.

Poster clamp strips are especially suitable for advertising presentations during fairs and in office and sales areas. The poster clamp strip is attached easily and rapidly to the wall or the ceiling. All that is needed is a nail and a piece of nylon string.

One packaging unit consists of 2 poster clamp strips made of aluminum.

**83891**Total width of the clamp strip:

610 mm



## UV-Polymerization / Venting



### Light Pistol Megalux CS

High-performance halogen UV-curing unit for the polymerization of UV-curing modeling pastes (e.g. Xpress) for shaping vents or executing repairs. Due to integrated program settings, curing time can be set to 20 or 40 seconds or continuous operation. Scope of delivery includes a small light guide. Its diameter of 3 mm is ideally suited for shaping vents. In addition, delivery contains the light safety shield and a spare halogen bulb 75 W.

**230 V**  
**115 V**

**51550**  
**51551**

*Output:* 130 W

*Spectral power distribution:* 400 – 500 nm

*Dimensions:* W 160 x H 80 x D 220 mm

*Spare parts:*

**Spare Light Guide**, 3 mm

**51552**

**Spare Halogen Bulb**, 75 W

**51503**



Also suitable for Demetron and Optilux.







## Silicone Technique – Cold Curing



### egger AQUA

egger AQUA is a floatable, elastic silicone for manufacturing splash and noise protection earmolds. The cold-curing material is extremely tensile and easy to use. egger AQUA is offered in different colors and can be used with the indirect method using a negative mold (duplicating material or plaster), as well as with the direct method.















egger AQUA, especially developed for the <b>direct method with a shorter processing and curing time.</b> Also suitable for the indirect method in the laboratory. <i>Packaging unit:</i> 8 double cartridges à 48 ml		Shore A (Final degree of hardness after 1 h)
		<b>35 shore</b>
	blue	REF 36952
	red	REF 36951
	red/blue – combination pack with 4 cartridges each	REF 36950
	yellow	REF 36957
	green	REF 36958
	orange	REF 36959

**egger flex/AB**

egger flex/AB is a high-quality and permanently elastic silicone for manufacturing tensile, soft-flexible BTE and noise protection earmolds. The different colors and shores are available and can be mixed with each other (e. g. for a soft auditory canal tip with a hard concha) egger flex/AB is usually processed with the indirect method using a negative mold (duplicating material or plaster). For a quick work process, we also offer the material with special formula for the direct use. Available colors: transparent and pink-transparent.






egger flex/AB materials, developed for the <b>indirect use in the laboratory with prolonged processing and curing times.</b> <i>Packaging unit:</i> 8 double cartridges à 48 ml		Shore A (Final degree of hardness after 1 hour)		
		<b>25 shore</b>	<b>40 shore</b>	<b>60 shore</b>
	reddish-transparent	REF 20801	REF 20810	REF 20741
	blue	REF 20806	REF 20704	REF 20743
	red	REF 20807	REF 20703	REF 20742
	colorless-transparent	REF 20735	REF 20818	REF 20740
	yellow		REF 20709	
	green		REF 20708	
	orange		REF 20719	
	beige		REF 20702	
	black		REF 20716	
	white		REF 20707	

egger flex/AB materials, especially developed for the <b>direct method with shorter process and curing times.</b> <b>Also suitable for the indirect method in the laboratory.</b> <i>Packaging unit:</i> 8 double cartridges à 48 ml		Shore A (Final degree of hardness after 1 hour)		
		<b>25 shore</b>	<b>40 shore</b>	<b>60 shore</b>
	reddish-transparent	REF 20701	REF 20700	REF 20705
	colorless-transparent		REF 20718	

**egger flex S-Pro – new formula**

egger flex S-Pro is a floatable, elastic silicone for splash and noise protection earmolds. The cold-curing material is very tensile and easy to use. Due to its extreme softness, the two components can be mixed easily. egger flex S-Pro is mixed manually and used in the egger impression syringe (REF 25500). The material can be used for the indirect method with a negative mold (duplicating material or plaster) or for the direct method.



egger flex S-Pro, especially developed for the <b>direct method with short processing and curing time.</b> Also suitable for the indirect method in the laboratory.		Shore A (Final degree of hardness after 1 h)	
		<b>40 shore</b>	
		<b>2 x 100 ml (Set incl. 1 pair of measuring spoons)</b>	<b>2 x 5 l</b>
	blue	REF 36985	REF 36990
	red	REF 36986	REF 36991
	beige	REF 36987	REF 36992

**Grip Stencil Set**

**36899**

Grip stencil for producing grip hollows on silicone earmolds, especially suited for hearing protection and splash protection earmolds. Can be used with the indirect and direct method for silicone cold-curing materials. Set includes a small (Ø 12 mm) and a large grip stencil (Ø 15 mm).



## Acrylic Technique



### Polymat, 230 V

A pressure polymerization unit for the professional laboratory! Due to the amplified temperature range up to 120 °C (248 °F), cold-curing and hot-curing acrylics can be processed quickly. Also suitable for the manufacturing of bubble-free negative molds of silicone. The curing process of the material can be adjusted by stepless electronic temperature control as well as pressure (3–6 bar) and time (up to 60 min.) control. The outer heater is very durable and facilitates the cleaning of the corrosion-resistant stainless steel pressure pot. Water is withdrawn by a drain. A high safety level is provided by a pressure controller, pressure relief valve, pressure switch, temperature limiter and safety interlock. Electricity and compressed air supply required.

*Pot diameter:* 220 mm

*Dimensions:* W 380 x D 330 x H 330 mm

*Inner height of pot:* 190 mm

*Weight:* 12 kg

*Heating output:* 1500 W

54200



## Laboratory Equipment



### Ultrasonic Cleansing Kit with Cedis Profi UltraActiv

The ideal starter kit for every hearing aid professional! The kit includes a compact, mobile ultrasonic cleansing unit with a wireless tank with handle! As the unit does not require a water supply at its location, the cleaner offers you unlimited mobility and is therefore ideally suited for every sales counter, service desk or fitting room. The inner container of stainless steel has a capacity of 750 ml. 5 time periods can be selected in the digital display. Includes basket. Additionally, the set comes with the disinfection concentrate UltraActiv (1 bottle of 500 ml for re-ordering: REF 86208) of our Cedis Profi range. With this especially for ultrasonic cleansing devices developed concentrate with antibacterial effectiveness you will achieve best results. For the exact dosing, the kit also contains one 100 ml dosing beaker and two 50 ml glass beakers for the ultrasonic bath.

86360 **No longer available.**

**Order**

Fast – easy – reliable.

The egger delivery service.

You order.

We deliver.

In time.


To the desired address.

“egger Labortechnik – How can I help you?”

You decide how you want to place your order.

The complete laboratory technique program including care products.

5 ways – one solution



Phone international: +49 831 5811360

Phone Hearing Protection: +49 831 5811390

Phone Earmold Production: +49 831 5811340

Monday to Friday 7:30 a.m. – 18:00 p.m.



Fax international: +49 831 5811314

Fax Earmold Production: +49 831 5811349



E-Mail international: [sales@egger-labor.de](mailto:sales@egger-labor.de)

E-Mail Earmold Production: [otoplastik@egger-labor.de](mailto:otoplastik@egger-labor.de)



egger Otoplastik + Labortechnik GmbH

Aybühlweg 59

87439 Kempten/Germany



Online-Shop

[www.egger-labor.com](http://www.egger-labor.com)

Delivery is effected according to our general terms and conditions. Technical modifications reserved.



